

CF3, the forward sequencing primer

Catalog No.	Size
GE100008	100 pmoles

Description

CF3, the forward sequencing primer for pCas-Guide system to sequence target sequence cloned in this vector.

Primer sequence: 5'-ACGATACAAGGCTGTTAGAGAG-3'

Appearance

Lyophilized. Please resuspend in 10 μ L H₂O to make a 10 μ M solution.

Storage

Store at 4°C.

Stability

Stable for a minimum of one year from date of receipt when stored at 4°C.

Reconstituted solution is stable for up to 4 weeks when stored at 4°C.

Related Products

pCas-Guide cloning kit

pCas-Guide plasmid

pCas-Scramble

Oligo annealing buffer

Complete genome editing tools OriGene offered, www.origene.com/cas9

Catalog No.

GE100001

GE100002

GE100003

GE100007

For Research Use Only. Not for use in diagnostic procedures.