



Transfecting immune cells with Viromers

Primary human macrophages transfection




primary M-2 macrophages **siRNA transfection** with **Viromer® BLUE** (right) or **GREEN** (left)






Human monocyte derived macrophages **plasmid transfection** with **Viromer Red**

Primary human monocytes ---- mRNA transfection



Human primary monocyte **GFP mRNA transfection** with **Viromer Red**

RAW 264.7 (mouse macrophages) transfection




siRNA transfection with Viromer[®] BLUE

- target gene AHA-1
- filled symbols : AHA-1 siRNA
- open symbols : control siRNA

GFP plasmid transfection with Viromer Red

THP-1 (monocytes) transfection



siRNA transfection using Viromer Blue (left) or Viromer Green (right)

Viromer products used



Cat# [TT100300](#)



Cat# [TT100301](#)



Cat# [TT100302](#)